


ANIMAL OPERATING TABLE

Features:-

1. Hydraulic up and down function.
2. IT is suitable for small animal's operation, such as sheep, dog etc.
3. It is reasonable structured, tight, high performance, easy operation and has very strong practicability.
4. Fully made of Premium Grade 304 Stainless Steel, it is heat-resistant, corrosion protective and Rust proof.
5. The surface of the table-top is V shaped and is smooth and easy to Clean and disinfect, thus, it can avoid cross infection.
6. Equipped with fluid collecting through system.



MODEL: - MMD-HVOT

Technical Data:-

TABLE TOP SIZE	24"× 62"
SUITABLE FOR ANIMALWEIGHT	UPTO 125KG
TABLE HEIGHT ADJUSTMENT BY HYDRAULIC [PUMP]	31"-40"
FLEX	90°
REFLEX	140°
SIDE TILT(BOTH SIDE)	20-25°
TENDLEBURG & REVERSETENDLEBURG	30°

MASPPO MEDICAL DEVICES

B-67/12, MAYAPURI INDL. AREA PHASE – I, NEW DELHI – 64

Mob. – 09871964472, 09958338559, E-Mail - masppo_medicaldevices@yahoo.com.

Web site:-www.masppomedicaldevices.com