

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:

MODEL NO: MMD3

Type Of Ceiling OT Light	Single Dome
Intensity	110000 Lux \pm 10%
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	57
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:

MODEL NO: MMD4

Type Of Ceiling OT Light	Single Dome
Intensity	130000 Lux \pm 10%
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	76
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:


MODEL NO: MMD5

Type Of Ceiling OT Light	Single Dome
Intensity	160000 Lux \pm 10%
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	95
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:

MODEL NO: MMD 3+3

Type Of Ceiling OT Light	Double Dome
Intensity	110000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	57 Each
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:


MODEL NO: MMD 4+4

Type Of Ceiling OT Light	Double Dome
Intensity	130000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	76 Each
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



Technical Data:

MODEL NO: MMD 5+5

Type Of Ceiling OT Light	Double Dome
Intensity	160000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	95 Each
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

With high light intensity and large spot size ,1200 mm. Depth of field provides good lighting without having to refocus when the surgical wound location changes. The light head can be positioned for away from the wound to ensure better laminar flow.



MODEL NO: MMD 48+48

Technical Data:

Type Of Ceiling OT Light	Double Dome
Intensity	160000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Colour Temperature	4000-5000k
Colour Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	48 Each
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	Aluminum
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional

Ceiling OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

With high light intensity and large spot size ,1200 mm. Depth of field provides good lighting without having to refocus when the surgical wound location changes. The light head can be positioned for away from the wound to ensure better laminar flow.



MODEL NO: MMD 48

Technical Data:

Type Of Ceiling OT Light	Single Dome
Intensity	160000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	48
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	Aluminum
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional