

Portable OT Light

Special features

High quality stenzilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

Technical Data:

Type Of Ceiling OT Light	Single Dome
Intensity	110000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	57
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional



MODEL NO: MMD 3


Portable OT Light

Special features

High quality sterilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

Technical Data:

Type Of Ceiling OT Light	Single Dome
Intensity	130000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	76
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional



MODEL NO: MMD 4

Portable OT Light

Special features

High quality sterilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

Technical Data:

Type Of Ceiling OT Light	Single Dome
Intensity	160000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	95
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	ABS
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional



MODEL NO: MMD 5

Portable OT Light

Special features

High quality sterilizable handle, with protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

Technical Data:

Type Of Ceiling OT Light	Single Dome
Intensity	160000 Lux $\pm 10\% * 2$
Size of Light Field	12-30cm
Color Temperature	4000-5000k
Color Reduction Index	95 RA
LED Life	<50000 Hrs.
Number Of LED	48
Diameter Of Light	500mm
Brightness Control	Capacitive Touch Panel
Material Of Head Lamp	Aluminum
Power Supply	220V/ 50HZ AC
Focusing	Adjustable
Battery Backup	Optional



MODEL NO: MMD 48